


VERTICAL ELECTROPHORESIS SYSTEM MP-8001



FEATURES

- SDS-PAGE gel setup
- Easy to use electrode assembly and leak-free.
- Four gel capacity
- Each module fits in to same buffer tank
- Run 4 gels in less than an hour
- Clearly display the running status
- Design of casting stand allows access of two gels simultaneously
- Casting frame with clamps
- Silicon rubber foaming against leakage
- Plastic combs
- Glass plate, spacer (thick) and short (thin) plate

TECHNICAL SPECIFICATION	
Gels	Two or Four
Glass plate size WxL	Thin plate 10.1x7.3cm and thick plate 10.1x8.2cm
Gel size WxL	Handcast 8.3x7.3cm and Precast 8.5x6.8cm
Buffer volume for two gel	650-750ml
Buffer volume for four gel	1000ml
Run time for SDS-PAGE	40-50min at 200V constant
Dimensions	12x16x18cm
Weight	1.0Kg