

WESTERN TRANSFER SYSTEM

ALS-WBT002



FEATURES

- Easily transferred SDS-PAGE gel with sample to the membrane
- Perform low intensity transfer by running at low voltage overnight
- Cooling unit prevents buffer depletion during high-intensity transfer
- Two gel capacity
- Design protects platinum wire electrode

TECHNICAL SPECIFICATION	
Cassettes	2
Transfer type	Tank
Blotting area	7.5x10cm
Fresh Buffer	450ml
Electrode distance	4cm
Transfer time	60min
Cooling	Blue ice cooling unit