

ELISA READER

ALS-INS007



KEY FEATURES

- Large LCD display, menu designed for convenient use.
 - Whole plate washing or single strip washing.
 - 12-way and 8-way manifold included.
 - Low residual volume by two pipettes.
 - Complete bottom washing.
 - Shaking and dipping function.
 - Automatic monitoring of vacuum and pressure, automatic rinse cycle.
 - Flat, V-bottom, or U-bottom plate and strip washing.
 - Large memory to store up to 48 user-programmed wash protocols.
- Emergency stop available and liquid warning.

TECHNICAL SPECIFICATION	
Manifold	8 pins and 12 pins
Wash time	0-99 times adjustable
Wash strips	1-12 adjustable
Dipping time	0-3600s adjustable
Resolution	50-3000ul/hole
Dispense precision	$C_v \leq 2\%$
Sipping time	0.1-9.9s
Power Supply	198-242V, 49-51Hz