

LIQUID NITROGEN CONTAINER

ALS-INS022



TECHNICAL SPECIFICATION

Specification	Unit	I	II	III	IV	V
Capacity	L	65	65	115	140	175
Empty weight	kg	30	30	40	45	58
Neck Diameter	mm	216	216	216	216	216
Outer Diameter	mm	681	681	681	681	681
Height	mm	700	726	796	910	1026
Static Evaporation L M	L/d	0.78	0.78	0.94	0.96	0.95
Static Holding Time	d	83	98	122	146	184
Number of Rocks	ea	6	6	6	6	6
Measurement of Rocks	mm	142×144	142×144	142×144	142×144	142×144
Measurement of box	mm	134×134	134×134	134×134	134×134	134×134
No. of stories for rocks	ea	4	5	6	8	10
10*10 cells for story	ea	100	100	100	100	100
Number of vials	ea	2400	3000	3600	4800	6000