

PH METER

ALS-INS031



FEATURES

- LCD displaying
- Manual/auto temperature compensation
- Function three-point calibration attached with pH composite electrode.

TECHNICAL SPECIFICATION		
Model		PH-METER31
Parameters Measured		pH/mV(ORP)/°C
Range	pH	(0.00~14.00) pH
	mV	(-1599.9)mV
	Temp	0.0~60.0°C
Resolution	pH	0.01pH
	mV	1mV
	Temp	0.1°C
Accuracy	pH	±0.02pH±1bit
	mV	±1%(FS)
	Temp	±0.1°C±1bit
Input resistance		>1×10 ⁹ Ω
Input Current		>1×10 ⁻¹⁰ A
Stability		(±0.05pH±1bit)/3h