

VISCOSITY METER

ALS-INS051



DESCRIPTION

A viscometer, viscosity meter or rheometer is a measuring instrument used to determine a fluid's internal flow resistance or viscosity.

KEY FEATURES

- Direct reading on display
- Big Ultra Bright Backlight LCD Display
- Auto Range function
- Temperature display
- Printer & Computer Interface (Model 200)
- Data Logger Software (Optional)not applicable in Model 60 & 100
- Display Shows selectable speed and Torque %

TECHNICAL SPECIFICATIONS

TECHNICAL SPECIFICATIONS					
Model No.	ALS-INS051-I	ALS-INS051-II	ALS-INS051-III	ALS-INS051-IV	ALS-INS051-V
Measuring Range mpa.s	20-2,000,000	20-2,000,000	80-40,000,000	20-6,000,000	40-40,000,000
No of Standard Spindle	4 0.3 to 60 8	4	6	4	6
Speed		0.3 to 100		0.1 to 200	
No of Speed		18		54	
Temperature Monitoring	Yes + / - 1%				
Measurement Accuracy	(FS) + / - 0.5 %				
Repeatability Power	(FS) 220 Ac /				
Supply	50Hz				